INSTRUCTIONS FOR USE

Route 92 Medical® Tenzing® 7 Delivery Catheter

Manufacturer

Route 92 Medical, Inc.
155 Bovet Road
Suite 100
San Mateo CA USA
customerservice@r92m.com


Emergo Europe B.V.
Prinsessegracht 20
2514 AP The Hague
The Netherlands

Refer to www.r92m.com/IFU for Summary of Safety & Clinical Performance or electronic instructions for use.

Device Description

The Route 92 Medical Tenzing 7 Delivery Catheter is a single-lumen, variable stiffness catheter with a long, tapered tip delineated by two radiopaque markers. The Delivery Catheter has a hydrophilic coating to increase lubricity. The proximal end has a luer hub. The Delivery Catheter is designed specifically for use with compatible catheters.

Indications for Use (United States FDA)

The Route 92 Medical Tenzing 7 Delivery Catheter is indicated for use with compatible catheters in facilitating the insertion and guidance of catheters into a selected blood vessel in the neurovascular system.

Intended Purpose

The Tenzing 7 Delivery Catheter is intended to be used to deliver large-bore catheters with an inner diameter of 0.068" or greater to the neurovasculature for the intended population of adult males and females under fluoroscopy using standard endovascular techniques. If required, a standard 0.014" or 0.016" neurovascular guidewire may be inserted through the Tenzing 7 Delivery Catheter.

Complications

Procedures requiring percutaneous catheter introduction should only be performed by physicians familiar with possible complications. Possible complications include but are not limited to the following: allergic reaction and anaphylaxis from contrast media; acute occlusion; air embolism; arteriovenous fistula; death; device malfunction; distal embolization; emboli; false aneurysm formation; hematoma or hemorrhage at access site; infection; intracranial hemorrhage; ischemia; kidney damage from contrast media; neurological deficits including stroke; vasospasm; and vessel perforation or dissection.

Compatibility

The Route 92 Medical Delivery Catheter is intended to be used with catheters with an inner diameter of 0.068" or greater and 0.014" or 0.016" diameter guidewires.

Warnings

- Do not advance or retract catheter against resistance without careful assessment of cause using fluoroscopy. If cause cannot be determined, withdraw catheter. Movement against resistance may result in catheter damage or patient injury.
- Do not use a device that has been damaged in any way. Use of a damaged device may result in complications.
- The catheter should only be used by physicians trained in interventional neuro-endovascular techniques.
- Testing has been limited to contrast media and saline. The use of this catheter for delivery of other solutions is not recommended.
- Do not use catheter with stents, retrievers, occlusion coils, glue, glue mixture or non-adhesive liquid embolic agent.

Angiography and Fluoroscopy Precautions

X-ray exposure from fluoroscopy and angiography poses risks including alopecia, burns ranging from skin reddening to ulcers, cataracts and delayed neoplasia. The probability of these risks occurring increases as procedure time and number of procedures increase. Care should be taken to minimize the X-ray radiation exposure of the patient and the operator by using sufficient shielding, reducing fluoroscopy times and modifying X-ray imaging techniques whenever possible.

Precautions

- Do not use high-powered contrast injection equipment. Use could result in damage to the device or vessel.
- Ensure target vessel diameter is appropriate and can accommodate catheter.
- Do not reuse or resterilize. The device is intended for single use only. Structural integrity and/or function may be impaired through reuse or cleaning.
- Store in cool, dry, dark place.
- Do not use opened or damaged packages.
- Use prior to the "Use By" date.
- Upon removal from package, inspect device to ensure it is not damaged.
- Do not expose device to solvents.
- Use device with fluoroscopic visualization and proper anti-coagulation agents.

- Hydrate catheter with heparinized saline before use. Once hydrated, do not allow the catheter to dry.
- Torquing the catheter while kinked may cause damage which could result in separation of the catheter shaft.
- Maintain a constant infusion of appropriate flush solution.
- If intraluminal device becomes lodged in catheter, or if the catheter becomes severely kinked, withdraw the entire system (intraluminal device, catheter and introducer sheath).
- Use only a steam source to shape the catheter tip.
- After steam shaping, inspect the catheter tip for damage. Do not use a catheter that has been damaged.
- To avoid damaging the catheter tip, do not steam shape the catheter tip more than twice.

Directions for Use

- The catheter tip may be steam shaped as follows:
 - Bend the tip to the preferred shape and hold above a steam source for 10 seconds or less.
 - Allow the tip to cool and inspect the tip for damage. Do not use if tip is damaged.
- Flush the Delivery Catheter hoop with heparinized saline to activate the hydrophilic coating. Flush inner diameter with saline.
- Carefully remove catheter from package hoop and inspect thoroughly, assuring that the device is not kinked or otherwise damaged. If there is any damage, replace with a new device. NOTE: Once hydrated, do not allow catheter to dry. Place catheter in basin of heparinized saline to keep hydrated.
- Attach an RHV to the Delivery Catheter. As needed insert a guidewire through the Delivery Catheter until the guidewire tip is just proximal to the tip of the Delivery Catheter.
- Attach a heparinized saline infusion line to the Delivery Catheter RHV.
- Insert the Delivery Catheter through the compatible catheter so that the tapered section of the Delivery Catheter extends past the tip of the catheter.
- Use standard catheterization techniques to advance the compatible catheter, the Delivery Catheter and guidewire to the targeted location under fluoroscopy. Remove the guidewire and Delivery Catheter.

Symbols

| | |
|--|--|
| | Caution |
| | Do not reuse |
| | Do not resterilize |
| | Consult instructions for use |
| | Lot number |
| | Catalogue number |
| | Sterilized using ethylene oxide |
| | Non-pyrogenic |
| | Use by date |
| | Manufacturer |
| | Prescription only - device restricted to use by or on the order of a physician |
| | Do not use if package is damaged |
| | Not made with natural rubber latex |
| | Keep dry |
| | Avoid temperature extremes |
| | Keep away from sunlight |
| | Medical Device |
| | Unique Device Identification |
| | Single Sterile Barrier System with Protective Packaging Inside |
| | Importer |

LIMITED WARRANTY

Route 92 warrants only that, at the time of manufacture, the Route 92 Medical Delivery Catheter (the "Products"): (i) are manufactured in accordance with good manufacturing practices, as required by the United States Food and Drug Administration; (ii) are true to label; (iii) conform to specifications; (iv) are free from defects in materials and workmanship; and (v) are not adulterated or misbranded within the meaning of the Federal, Food, Drug and Cosmetic Act as amended. Route 92 does not warrant a good effect or against any ill effect following the Products' use, and Route 92 makes no warranties that the Products are effective under all circumstances. THE FOREGOING WARRANTIES ARE EXCLUSIVE,

INSTRUCTIONS FOR USE

Route 92 Medical® Tenzing® 7 Delivery Catheter

ARE IN LIEU OF, AND SUPERSEDE ALL OTHER WARRANTIES OF ANY KIND, WHETHER WRITTEN, ORAL, EXPRESS, OR IMPLIED. THERE ARE NO WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR NONINFRINGEMENT. No representative of Route 92 may change any of the foregoing. Buyer's sole and exclusive remedy for any claim arising by reason of or in connection with the sale, purchase, delivery or use of Products, regardless of whether such claim is based on tort law, breach of contract, breach of warranty or any other legal theory shall be, at Route 92's option, the repair or replacement of any nonconforming Products.

Serious Incidents

In the event a serious incident has occurred in relation to the use of this device notify customerservice@r92m.com and applicable competent authorities.

Patents: www.r92m.com/patents

Symbols Glossary: www.r92m.com/symbols

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