INSTRUCTIONS FOR USE

Route 92 Medical® HiPoint™ 70 and HiPoint 88 Reperfusion System

Manufacturer
Route 92 Medical, Inc.
155 Bovet Road
Suite 100
San Mateo, CA USA
customerservice@r92m.com

Emergo Europe B.V.
Westervoortse dijk 60
6827 AT Arnhem
The Netherlands

Refer to www.r92m.com/IFU for Summary of Safety & Clinical Performance and electronic instructions for use.

Device Description

A Delivery Catheter and Aspiration Catheter comprise each Route 92 Medical Reperfusion System. The Delivery Catheter is a single-lumen, variable stiffness catheter with a long, tapered tip delineated by two radiopaque markers. The proximal end has a luer hub. The Delivery Catheter is designed specifically for use with the Aspiration Catheter. The Aspiration Catheter is a single-lumen, variable stiffness catheter with a radiopaque marker at the distal tip. Both the Delivery Catheter and the Aspiration Catheter are coated with a hydrophilic coating to facilitate movement. A Clip is provided to temporarily secure the Delivery Catheter to the Aspiration Catheter during delivery.

Figure 1 shows the Delivery Catheter inserted into the Aspiration Catheter.

Figure 1: Delivery Catheter inserted through Aspiration Catheter

Figure 2 shows the Aspiration Catheter inserted into a long sheath with the Delivery Catheter removed.

Figure 2: Aspiration Catheter inserted into Sheath

Intended Purpose

The Route 92 Medical HiPoint Reperfusion System is intended for use in the introduction of interventional devices into the neurovasculature and for aspiration of thrombus in ischemic stroke patients

Complications

Procedures requiring percutaneous catheter introduction should only be performed by physicians familiar with possible complications. Possible complications include but are not limited to the following: allergic reaction and anaphylaxis from contrast media; acute occlusion; additional surgical intervention: air embolism; arteriovenous fistula; death; device malfunction; distal embolization; emboli; false aneurysm formation; hematoma or hemorrhage at access site; sterile inflammation or granulomas at the access site; infection; intracranial hemorrhage; ischemia; kidney damage from contrast media; neurological deficits including stroke; residual thrombus; tissue necrosis, transient or long lasting: vasospasm; and vessel perforation or dissection.

Compatibility

- The HiPoint 70 Reperfusion System is compatible with catheters or sheaths with an inner diameter of 0.088" (2.24 mm).
- The HiPoint 88 Reperfusion System is compatible with sheaths with an inner diameter of 0.106" (2.69 mm).

Only guidewires may be introduced through the Route 92 Medical Delivery Catheters. The Delivery Catheters are not compatible with embolic coils, stent retrievers or other interventional devices.

The Route 92 Medical Delivery Catheters are compatible with guidewires 0.016" or less in diameter.

Warnings

- Do not advance or retract catheter against resistance without careful assessment of cause using fluoroscopy. If cause cannot be determined, withdraw catheter. Movement against resistance may result in catheter damage or patient injury.
- Do not use a device that has been damaged in any way. Use of a damaged device may result in complications.
- The system should only be used by physicians trained in interventional neuro-endovascular techniques.

Precautions

- **Do not use high-powered contrast injection equipment. Use could result in damage to the device or vessel.**
- Ensure target vessel diameter is appropriate and can accommodate catheter.
- Do not reuse or resterilize. The device is intended for single use only. Structural integrity and/or function may be impaired through reuse or cleaning.
- Store in cool, dry, dark place.
- Do not use opened or damaged packages.
- Use prior to the "Use By" date.
- Upon removal from package, inspect each device to ensure no damage.
- Do not expose device to solvents.
- Use device with fluoroscopic visualization and proper anti-coagulation agents.
- Hydrate catheter with heparinized saline before use. Once hydrated, do not allow the catheter to dry.

- Torquing the catheter while kinked may cause damage which could result in separation of the catheter shaft.
- Maintain a constant, pressurized, Heparinized saline infusion on all devices.
- If intraluminal device becomes lodged in catheter, or if the catheter becomes severely kinked, withdraw the entire system (intraluminal device, catheter and introducer sheath).
- When aspirating, aspirate for the minimum time required to remove thrombus.
- Monitor intra-procedural blood loss and manage as appropriate.
- Use only a steam source to shape the Delivery Catheter tip.
- After steam shaping, inspect the Delivery Catheter tip for damage. Do not use a catheter that has been damaged.
- To avoid damaging the Delivery Catheter tip, do not steam shape the catheter tip more than twice.
- In between aspiration passes, withdraw the Aspiration Catheter(s) from the patient and clean any residual thrombus prior to reinsertion and subsequent contrast injection.


Directions for Use

1. The Delivery Catheter tip may be steam shaped as follows:
 - Bend the tip to the preferred shape and hold above a steam source for 10 seconds or less.
 - Allow the tip to cool and inspect the tip for damage. Do not use if tip is damaged.
2. Flush the Aspiration Catheter hoop with heparinized saline to flush the inner diameter of the catheter and activate the hydrophilic coating.
3. Flush the Delivery Catheter hoop with heparinized saline to activate the hydrophilic coating. Flush inner diameter with saline.
4. Carefully remove devices from packaging and inspect thoroughly, ensuring the devices are not kinked or otherwise damaged. If there is any damage, replace with a new device. NOTE: Once hydrated, do not allow catheters to dry. Between passes, place catheters in a basin of heparinized saline to keep hydrated.
5. Place the guide/sheath at an appropriate location proximal to the occlusion, according to the manufacturer's instructions for use.
6. Attach a large-bore rotating hemostasis valve (RHV) to the guide/sheath. Attach a pressurized, heparinized saline infusion line with a 3-way stopcock to the RHV.
7. Using the 3-way stopcock, manually flush the sheath with a syringe. Repeat manual flush through the 3-way stopcock any time backbleeding into the RHV is observed.
8. Attach an RHV to the Delivery Catheter. As preferred, insert a guidewire through the Delivery Catheter until the guidewire tip is just proximal to the tip of the Delivery Catheter.
9. Attach a pressurized, heparinized saline infusion line to the Delivery Catheter RHV.
10. Insert the Delivery Catheter through the Aspiration Catheter so that the tapered section of the Delivery Catheter extends past the tip of the Aspiration Catheter.
11. The Clip may be used to temporarily secure the Delivery Catheter to the Aspiration Catheter during delivery. Place the Aspiration Catheter ribbon into the Clip groove first, followed by the Delivery Catheter.
12. Gently insert Reperfusion System (Delivery Catheter and Aspiration Catheter) through the RHV of the guide/sheath.
13. Use conventional catheterization techniques to advance the Reperfusion System to the thrombus under fluoroscopy. Position the Aspiration Catheter at the thrombus. Remove the guidewire (if used) and Delivery Catheter and Clip. As needed, the HiPoint 70 Reperfusion System may be inserted through the 088 Aspiration Catheter in order to access more distal vessels.
14. Attach vacuum source to the RHV of guide/sheath. Firmly tighten the RHV around the Aspiration Catheter to prevent backflow and entrainment of air during application of vacuum. Apply vacuum to aspirate thrombus.
15. If there is not free flow through the Aspiration Catheter at the end of the aspiration, remove the Aspiration Catheter from the patient while maintaining vacuum. Use a syringe filled with heparinized saline to flush the Aspiration Catheter to remove any obstructing thrombus.
16. In between aspiration passes and prior to contrast injection, withdraw the Aspiration Catheter(s) from the patient and clean any residual or obstructing thrombus from the Aspiration Catheter(s) prior to re-insertion and subsequent contrast injection.
17. Repeat aspiration procedure as necessary.

Symbols (see specific package labeling for applicable symbols)

	Caution
	Do not reuse
	Do not resterilize
	Consult instructions for use
	Lot number
	Catalogue number
	Sterilized using ethylene oxide
	Non-pyrogenic
	Use by date
	Manufacturer
	Prescription only - device restricted to use by or on the order of a physician

INSTRUCTIONS FOR USE**Route 92 Medical® HiPoint™ 70 and HiPoint 88 Reperfusion System**

	Do not use if package is damaged
	Not made with natural rubber latex
	Keep dry
	Avoid temperature extremes
	Keep away from sunlight
	Medical Device
	Unique Device Identification
	Single sterile barrier system with protective packaging outside

LIMITED WARRANTY

Route 92 warrants only that, at the time of manufacture, the Route 92 Medical Reperfusion System (the "Products"): (i) are manufactured in accordance with good manufacturing practices, as required by the United States Food and Drug Administration; (ii) are true to label; (iii) conform to specifications; (iv) are free from defects in materials and workmanship; and (v) are not adulterated or misbranded within the meaning of the Federal Food, Drug, and Cosmetic Act as amended. Route 92 does not warrant a good effect or against any ill effect following the Products' use, and Route 92 makes no warranties that the Products are effective under all circumstances. THE FOREGOING WARRANTIES ARE EXCLUSIVE, ARE IN LIEU OF, AND SUPERSEDE ALL OTHER WARRANTIES OF ANY KIND, WHETHER WRITTEN, ORAL, EXPRESS, OR IMPLIED. THERE ARE NO WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR NONINFRINGEMENT. No representative of Route 92 may change any of the foregoing. Buyer's sole and exclusive remedy for any claim arising by reason of or in connection with the sale, purchase, delivery or use of Products, regardless of whether such claim is based on tort law, breach of contract, breach of warranty or any other legal theory shall be, at Route 92's option, the repair or replacement of any nonconforming Products.

Serious Incidents

In the event a serious incident has occurred in relation to the use of this device notify customerservice@r92m.com and applicable competent authorities.

Patents: www.r92m.com/patents

Symbols Glossary: www.r92m.com/symbols

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