INSTRUCTIONS FOR USE

Route 92 Medical® FreeClimb™ 70 and FreeClimb 88 Reperfusion Systems

Manufacturer
Route 92 Medical
155 Bovet Road
Suite 100
San Mateo CA USA
customerservice@r92m.cc

Emergo Europe B.V.
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands

Refer to www.r92m.com/IFU for Summary of Safety & Clinical Performance and electronic instructions for use.

Device Description

A Delivery Catheter and Aspiration Catheter comprise the Route 92 Medical FreeClimb Reperfusion Systems. The Delivery Catheter is a single-lumen, variable stiffness catheter with a long, tapered tip delineated by two radiopaque markers. The proximal end has a luer hub. The Delivery Catheter is designed specifically for use with the Aspiration Catheter. The Aspiration Catheter is a single-lumen, variable stiffness catheter with a radiopaque marker at the distal tip. Both the Delivery Catheter and the Aspiration Catheter are coated with a hydrophilic coating to facilitate movement.

- The Route 92 Medical FreeClimb 88 Reperfusion System is provided in two lengths, 132 cm and 125 cm. The provided Delivery Catheter is the same for both lengths.
- The Route 92 Medical FreeClimb 70 Reperfusion System is provided in one length, 132 cm.

Intended Purpose

The Route 92 Medical FreeClimb Reperfusion Systems are intended for use in the introduction of interventional devices into the neurovasculature and for aspiration of thrombus in ischemic stroke patients.

Complications

Procedures requiring percutaneous catheter introduction should only be performed by physicians familiar with possible complications. Possible complications include but are not limited to the following: allergic reaction and anaphylaxis from contrast media; acute occlusion; additional surgical intervention: air embolism; arteriovenous fistula; death; device malfunction; distal embolization; emboli; false aneurysm formation; hematoma or hemorrhage at access site; sterile inflammation or granulomas at the access site; infection; intracranial hemorrhage; ischemia; kidney damage from contrast media; neurological deficits including stroke; residual thrombus; tissue necrosis, transient or long lasting; vasospasm; and vessel perforation or dissection.

Compatibility

The Route 92 Medical FreeClimb 88 Reperfusion System is compatible with catheters or sheaths with a minimum inner diameter of 0.106" (2.69 mm).

The Route 92 Medical FreeClimb 70 Reperfusion System is compatible with catheters or sheaths with a minimum inner diameter of 0.088" (2.24 mm).

Only guidewires may be introduced through the Route 92 Medical Delivery Catheter. The Delivery Catheter is not compatible with embolic coils, stent retrievers or other interventional devices.

The Route 92 Medical Delivery Catheters are compatible with guidewires 0.016" or less in diameter.

Warnings

- Do not advance or retract catheter against resistance without careful assessment of cause using fluoroscopy. If cause cannot be determined, withdraw catheter. Movement against resistance may result in catheter damage or patient injury.
- Do not use a device that has been damaged in any way. Use of a damaged device may result in complications.
- The system should only be used by physicians trained in interventional neuro-endovascular techniques.

Precautions

- **Do not use high-powered contrast injection equipment. Use could result in damage to the device or vessel.**
- Ensure target vessel diameter is appropriate and can accommodate catheter. Use device in vessels that are larger than the outer diameter of the device
 - FreeClimb 88 Aspiration Catheter: ≥2.7mm
 - FreeClimb 70 Aspiration Catheter: ≥2.2 mm
- Do not reuse or resterilize. The device is intended for single use only. Structural integrity and/or function may be impaired through reuse or cleaning.
- Store in cool, dry, dark place.
- Do not use opened or damaged packages.
- Use prior to the "Use By" date.
- Upon removal from package, inspect each device to ensure no damage.
- Do not expose device to solvents.
- Use device with fluoroscopic visualization and proper anti-coagulation agents.
- Hydrate catheter with heparinized saline before use. Once hydrated, do not allow the catheter to dry.
- Torquing the catheter while kinked may cause damage which could result in separation of the catheter shaft.
- Maintain a constant, pressurized, heparinized saline infusion on all devices.
- If intraluminal device becomes lodged in catheter, or if the catheter becomes severely kinked, withdraw the entire system (intraluminal device, catheter and introducer sheath).
- When aspirating, aspirate for the minimum time required to remove thrombus.
- Monitor intra-procedural blood loss and manage as appropriate.
- Use only a steam source to shape the Delivery Catheter tip.
- After steam shaping, inspect the Delivery Catheter tip for damage. Do not use a catheter that has been damaged.
- To avoid damaging the Delivery Catheter tip, do not steam shape the catheter tip more than twice.

- In between aspiration passes, withdraw the Aspiration Catheter from the patient and clean any residual thrombus prior to reinsertion and subsequent contrast injection.

Directions for Use


1. The Delivery Catheter tip may be steam shaped as follows:
 - Bend the tip to the preferred shape and hold above a steam source for 10 seconds or less.
 - Allow the tip to cool and inspect the tip for damage. Do not use if tip is damaged.
2. Flush the Aspiration Catheter hoop with heparinized saline to activate the hydrophilic coating. Flush inner diameter with heparinized saline.
3. Flush the Delivery Catheter hoop with heparinized saline to activate the hydrophilic coating. Flush inner diameter with heparinized saline.
4. Carefully remove devices from packaging and inspect thoroughly, ensuring the devices are not kinked or otherwise damaged. If there is any damage, replace with a new device. NOTE: Once hydrated, do not allow catheters to dry. Between passes, place catheters in a basin of heparinized saline to keep hydrated.
5. Place the guide/sheath at an appropriate location proximal to the occlusion, according to the manufacturer's instructions for use.
6. Attach a large-bore rotating hemostasis valve (RHV) to the guide/sheath. Attach a pressurized, heparinized saline infusion line with a 3-way stopcock to the RHV.
7. Using the 3-way stopcock, manually flush the sheath with a syringe. Repeat manual flush through the 3-way stopcock any time backbleeding into the RHV is observed.
8. Attach an RHV to the Delivery Catheter. As preferred, insert a guidewire through the Delivery Catheter until the guidewire tip is just proximal to the tip of the Delivery Catheter.
9. Attach a pressurized, heparinized saline infusion line to the Delivery Catheter RHV.
10. Attach an RHV to the Aspiration Catheter and attach a pressurized, heparinized saline infusion line to the Aspiration Catheter RHV.
11. Insert the Delivery Catheter through the Aspiration Catheter so that all of the tapered section of the Delivery Catheter extends past the tip of the Aspiration Catheter.
12. Gently insert Reperfusion System (Delivery Catheter and Aspiration Catheter) through the RHV of the guide/sheath.
13. Use conventional catheterization techniques to advance the Reperfusion System to the thrombus under fluoroscopy. Position the Aspiration Catheter at the thrombus. Remove the guidewire (if used), and Delivery Catheter.
14. Attach vacuum source to the RHV of guide/sheath. Firmly tighten the RHV around the Aspiration Catheter to prevent backflow and entrainment of air during application of vacuum. Apply vacuum to aspirate thrombus.
15. If there is not free flow through the Aspiration Catheter at the end of the aspiration, remove the Aspiration Catheter from the patient while maintaining vacuum. Use a syringe filled with heparinized saline to flush the Aspiration Catheter to remove any obstructing thrombus.
16. In between aspiration passes and prior to contrast injection, withdraw the Aspiration Catheter from the patient and clean any residual or obstructing thrombus from the Aspiration Catheter prior to re-insertion and subsequent contrast injection.
17. Repeat aspiration procedure as necessary.

Symbols (see specific package labeling for applicable symbols)

| | |
|--|--|
| | Caution |
| | Do not reuse |
| | Do not resterilize |
| | Consult instructions for use |
| | Lot number |
| | Catalogue number |
| | Sterilized using ethylene oxide |
| | Non-pyrogenic |
| | Use by date |
| | Manufacturer |
| | Prescription only - device restricted to use by or on the order of a physician |
| | Do not use if package is damaged |
| | Not made with natural rubber latex |
| | Keep dry |
| | Avoid temperature extremes |
| | Keep away from sunlight |
| | Medical Device |

INSTRUCTIONS FOR USE

Route 92 Medical® FreeClimb™ 70 and FreeClimb 88 Reperfusion Systems

| | |
|--|---|
|  | Unique Device Identification |
|  | Single sterile barrier system with protective packaging outside |

LIMITED WARRANTY

Route 92 warrants only that, at the time of manufacture, the Route 92 Medical Reperfusion Systems (the "Products"): (i) are manufactured in accordance with good manufacturing practices, as required by the United States Food and Drug Administration; (ii) are true to label; (iii) conform to specifications; (iv) are free from defects in materials and workmanship; and (v) are not adulterated or misbranded within the meaning of the Federal, Food, Drug and Cosmetic Act as amended. Route 92 does not warrant a good effect or against any ill effect following the Products' use, and Route 92 makes no warranties that the Products are effective under all circumstances. THE FOREGOING WARRANTIES ARE EXCLUSIVE, ARE IN LIEU OF, AND SUPERSEDE ALL OTHER WARRANTIES OF ANY KIND, WHETHER WRITTEN, ORAL, EXPRESS, OR IMPLIED. THERE ARE NO WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR NONINFRINGEMENT. No representative of Route 92 may change any of the foregoing. Buyer's sole and exclusive remedy for any claim arising by reason of or in connection with the sale, purchase, delivery or use of Products, regardless of whether such claim is based on tort law, breach of contract, breach of warranty or any other legal theory shall be, at Route 92's option, the repair or replacement of any nonconforming Products.

Serious Incidents

In the event a serious incident has occurred in relation to the use of this device notify customerservice@r92m.com and applicable competent authorities.

Patents: www.r92m.com/patents

Symbols Glossary: www.r92m.com/symbols

Fabricant

Route 92 Medical
155 Bovet Road
Suite 100
San Mateo, CA USA
customerservice@r92m.com

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