INSTRUCTIONS FOR USE

Route 92 Medical® Base Camp® Sheath System 2.0

Manufacturer
Route 92 Medical
155 Bovet Road
Suite 100
San Mateo, CA 94402 USA
info@r92m.com

Refer to www.r92m.com/IFU for Summary of Safety & Clinical Performance or electronic instructions for use.

Device Description

The Route 92 Medical Base Camp® Sheath System 2.0 is comprised of a Sheath, a Dilator, a Navigating Catheter and an RHV (rotating hemostasis valve). The Base Camp Sheath is a single-lumen, variable stiffness catheter with a radiopaque marker on the distal end. The inner lumen of the catheter is compatible with 8F or smaller catheters. The Dilator may be placed within the Sheath to facilitate percutaneous introduction of the Sheath into a femoral artery. The Dilator has a radiopaque marker at the distal tip. The Navigating Catheter is a single-lumen, variable stiffness catheter with a radiopaque marker at the distal tip. The Navigating Catheter is compatible with the Sheath and has a shaped distal end to facilitate placement. All of the catheters are coated with hydrophilic coating.

Hydrophilic Coating Length	
Sheath	8 cm from the distal end
Dilator	15 cm from the distal end
Navigating Catheter	40 cm from the distal end

Intended Purpose

The Route 92 Medical Sheath System 2.0 is indicated for the introduction of interventional devices into the neuro vasculature.

Contraindications

None known.

Complications

Procedures requiring percutaneous catheter introduction should only be performed by physicians familiar with possible complications. Possible complications include but are not limited to the following: allergic reaction and anaphylaxis from contrast media; acute occlusion; air embolism; arteriovenous fistula; death; device malfunction; distal embolization; emboli; false aneurysm formation; hematoma or hemorrhage at access site; sterile inflammation or granulomas at the access site; additional surgical intervention; tissue necrosis (transient or long-lasting); infection; intracranial hemorrhage; ischemia; kidney damage from contrast media; neurological deficits including stroke; thrombosis; vasospasm; and vessel perforation or dissection.

Compatibility

The Route 92 Medical Sheath is intended to be used with 8F or smaller catheters. An 0.035" guidewire may be used through the Dilator and the Navigating Catheter.

Dimensional Specifications			
Device Component	Working Length	Inner Diameter	Outer Diameter
Sheath	90 cm (35.4")	2.7 mm (0.106")	3.1 mm (0.122")
	80 cm (31.5")		
Dilator	103 cm (40.6")	1.0 mm (0.040")	2.6 mm (0.102")
Navigating Catheter	113 cm (44.5")	1.0 mm (0.040")	2.7 mm (0.104")

Warnings

- Do not advance or retract catheters against resistance without careful assessment of cause using fluoroscopy. If cause cannot be determined, withdraw catheter. Movement against resistance may result in catheter damage or patient injury.
- Do not use a device that has been damaged in any way. Use of a damaged device may result in complications.
- The catheters should only be used by physicians trained in interventional neuro-endovascular techniques.
- Testing has been limited to contrast media and saline. The use of this catheter for delivery of other solutions is not recommended.
- Do not use device in vessels that are smaller than the outer diameter of the device or within tortuous vasculature.
- Do not use high-powered contrast injection equipment. Use could result in damage to the device or vessel.
- Do not reuse or resterilize. The device is intended for single use only. Structural integrity and/or function may be impaired through reuse or cleaning.
- Do not use opened or damaged packages.
- Do not expose device to solvents.
- The Sheath System is coated with a hydrophilic coating at the distal end of each device for a length of 8 cm for the Sheath, 15 cm for the Dilator and 40 cm for the Navigating Catheter. Please refer to the Directions for Use section for further information on how to prepare and use the system to ensure it performs as intended. Failure to abide by the warnings in this labeling might result in damage to the device coating, which may necessitate intervention or result in serious adverse events.
- Use caution when manipulating, advancing and/or withdrawing these devices through needles, metal cannulas, stents, or other devices with sharp edges, or through tortuous or calcified blood vessels. Manipulation, advancement, and/or withdrawal past sharp or beveled edges may result in destruction and/or separation of the outer coating, which may lead to clinical adverse events, resulting in coating material remaining in the vasculature or device damage. This may result in adverse events requiring additional intervention.

Angiography and Fluoroscopy Precautions

X-ray exposure from fluoroscopy and angiography poses risks including alopecia, burns ranging from skin reddening to ulcers, cataracts and delayed neoplasia. The probability of these risks occurring increases as procedure time and number of procedures increase. Care should be taken to minimize the X-ray radiation exposure of the patient and the operator by using sufficient shielding, reducing fluoroscopy times and modifying X-ray imaging techniques whenever possible.

Precautions

- Store in cool, dry, dark place.
- Use prior to the "Use By" date.
- Upon removal from package, inspect device to ensure it is not damaged.
- Use device with fluoroscopic visualization and proper anti-coagulation agents.
- Hydrate each catheter with heparinized saline before use. Once hydrated, do not allow the catheter to dry.
- Torquing the catheter while kinked may cause damage which could result in separation of the catheter shaft.
- Maintain a constant infusion of appropriate flush solution.
- If intraluminal device becomes lodged in catheter, or if the catheter becomes severely kinked, withdraw the entire system (intraluminal device, catheter and introducer sheath).
- Avoid wiping the device with dry gauze as this may damage the device coating.
- Avoid excessive wiping of the coated device
- Avoid using alcohol, antiseptic solutions, or other solvents to pre-treat the device because this may cause unpredictable changes in the coating which could affect the device safety and performance.
- Avoid pre-soaking the devices as this may impact the coating performance.

Directions for Use

- Carefully remove the catheters from the packaging and inspect thoroughly, assuring that none of the devices are kinked or otherwise damaged. If there is any damage, replace with a new device.
- Flush the inner lumens of the catheters with saline.
- Fully insert the Dilator into the Sheath. The luer of the Dilator is designed to press fit into the Sheath luer.
- If not using an introducer sheath
 - Access artery using standard techniques
 - Advance the Sheath/Dilator assembly over a guidewire in the artery to the intended position and remove the dilator
- If using an introducer sheath
 - Place an appropriately sized introducer sheath in the artery using standard techniques
 - Advance the Sheath/Dilator assembly through the introducer sheath to the intended position and remove the Dilator
- Manually flush the Sheath and attach an RHV with a heparinized saline flush line to the Sheath.
- Insert the Navigating Catheter through the Sheath to assist in placement of the Sheath.
- Advance the Sheath to the intended location under fluoroscopic guidance using standard endovascular techniques.
- Place compatible catheters through the Sheath as needed.
- When use of the Sheath is complete, remove the Sheath using standard techniques.
- Dispose according to standard hospital practices.

Symbols

	Caution
	Do not reuse
	Do not resterilize
	Consult instructions for use
	Lot number
	Catalogue number
	Sterilized using ethylene oxide
	Non-pyrogenic
	Use by date
	Manufacturer
	Prescription only - device restricted to use by or on the order of a physician
	Do not use if package is damaged
	Not made with natural rubber latex
	Keep dry

	Keep away from sunlight
	Medical Device
	Unique Device Identification
	Single sterile barrier system with protective packaging outside
	Single sterile barrier system

LIMITED WARRANTY

Route 92 warrants only that, at the time of manufacture, the Route 92 Medical Base Camp Sheath System 2.0 (the "Products"): (i) are manufactured in accordance with good manufacturing practices, as required by the United States Food and Drug Administration; (ii) are true to label; (iii) conform to specifications; (iv) are free from defects in materials and workmanship; and (v) are not adulterated or misbranded within the meaning of the Federal, Food, Drug and Cosmetic Act as amended. Route 92 does not warrant a good effect or against any ill effect following the Products' use, and Route 92 makes no warranties that the Products are effective under all circumstances. THE FOREGOING WARRANTIES ARE EXCLUSIVE, ARE IN LIEU OF, AND SUPERSEDE ALL OTHER WARRANTIES OF ANY KIND, WHETHER WRITTEN, ORAL, EXPRESS, OR IMPLIED. THERE ARE NO WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR NONINFRINGEMENT. No representative of Route 92 may change any of the foregoing. Buyer's sole and exclusive remedy for any claim arising by reason of or in connection with the sale, purchase, delivery or use of Products, regardless of whether such claim is based on tort law, breach of contract, breach of warranty or any other legal theory shall be, at Route 92's option, the repair or replacement of any nonconforming Products.

Serious Incidents

In the event a serious incident has occurred in relation to the use of this device notify customerservice@r92m.com and applicable competent authorities.

Patents: www.r92m.com/patents
Symbols Glossary: www.r92m.com/symbols

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